

Mass propagation of *Nicotiana tabacum* cv. BY-2 suspension cells in the Finesse SmartRocker bioreactor

Nicole Imseng, Stephan C. Kaiser, Nadezda Perepelitsa, and Dieter Eibl

Zurich University of Applied Sciences, School of Life Science and Facility Management, Institute of Biotechnology, Biochemical Engineering and Cell Cultivation Techniques, Switzerland

Summary

This application note describes the mass propagation of *Nicotiana tabacum* cv. BY-2 (cultivar tobacco Bright Yellow-2) suspension cells in the Finesse SmartRocker bioreactor at laboratory scale. Maximum fresh cell weights of about $351 \pm 9 \text{ g} \cdot \text{L}^{-1}$, corresponding to packed cell volumes of $74 \pm 2 \%$, were achieved, which shows the applicability of the Finesse SmartRocker bioreactor to cultivate fast growing plant suspension cells.

Introduction

The fast-growing tobacco *Nicotiana tabacum* cv. BY-2 suspension cell line has been established as model cell line in many laboratories and is one of the most widely used plant cell cultures [Kirchhoff, 2012]. Successful cell expansion in wave-mixed bioreactors at 1 L scale have previously been demonstrated [Raven *et al.*, 2010; Brändli *et al.*, 2012; Werner *et al.*, 2010], where packed cell volumes (pcv) of up to 70 % have been achieved. Using initial cell concentrations of between 10 and 15 % pcv, a three-day exponential growth period with growth rates in the range of 0.54 and 0.79 d^{-1} has been observed.

In the present study the suitability of the Finesse SmartRocker bioreactor (see Fig. 1) for the cultivation of BY-2 plant suspension cells is demonstrated. The SmartBag bioprocess container with a maximum working volume of 5 L is equipped with single-use sensors for pH and DO measurement (Finesse TruFluor). They are operated via the G3Lab Universal controller, which incorporates SmartMFCs for up to six gases, a standard set of four peristaltic SmartPumps, a SmartMotor driver and SmartTransmitters for TruFluor or electrochemical sensors. The system is connected

to an industrial DeltaV enabled TruLogic™ controller operated with TruBio® DV software for bioreactor control, data acquisition and process visualization. The cells were cultivated in batch-mode using chemically defined MSMOplus (3% sucrose) culture medium.

In order to evaluate the reproducibility, three individual cultivations with each running 10 days were realized. Maximum fresh cell weights of approximately $350 \text{ g} \cdot \text{L}^{-1}$, corresponding to packed cell volumes of 74 %, were achieved. Even though the culture becomes increasingly viscous during the cultivation ($>100 \text{ mPa} \cdot \text{s}$, based on own studies), increasing the rocking rate was sufficient to guarantee well-mixed conditions and high cell growth.



Fig. 1: Picture of the Finesse SmartRocker bioreactor with PC based Finesse G3Lab Universal controller

Materials and methods

Cell line and culture media

For the experiments, the tobacco suspension cell line of *Nicotiana tabacum* cv. Bright Yellow-2 (BY-2) (established and kindly provided by the group of Dr. Stefan Schillberg, IME Aachen, Germany) was used. The cultivations were realized using MSMOplus medium (Cell Culture Technologies, CH). The composition is summarized in Tab. 1. The pH value was adjusted to pH 5.8 using 1 M sodium hydroxide solution.

Tab. 1: Composition of the MSMOplus culture medium.

Component	Concentration
MSMO macroelements	
MSMO microelements	
Inositol	100 mg·L ⁻¹
Thiamine HCl	1 mg·L ⁻¹
Sucrose	30 g·L ⁻¹
Monosodium phosphate	230 mg·L ⁻¹
Kanamycin	0.1 mg·L ⁻¹
2,4-Dichlorophenoxyacetic acid	0.5 mg·L ⁻¹

Overview of setting-up procedures

Day -7: Inoculum production with *Nicotiana tabacum* cv. BY-2 suspension cells in shake flasks.

Day 0: Medium filling and equilibration, sensor calibration. Inoculation and starting-up of the SmartBag bioprocess container with a seeding concentration of 10 % pcv in 5.0 L MSMOplus medium (Cell Culture Technologies, CH), Sampling.

Days 1 – 10: Sampling, analytics and harvest.

Inoculum preparation

The inoculum for the Finesse SmartRocker bioreactor was produced in five single-use shake flasks (Corning) at a maximum working volume of 300 mL. The cells were inoculated with cell concentrations of 15 % pcv. The flasks were shaken at a shaking rate of 180 rpm and amplitude of 25 mm (Multitron shaker, Infors HT, CH) in darkness. After three days, the cell suspension had a pcv of 50 % and cell viability was assessed to be above 95 %.

Bioreactor preparation and inoculation

About one hour prior to inoculation, a Flexboy® media bag (Sartorius Stedim Biotech) with 3 L MSMOplus culture medium was connected to the SmartBag bioprocess container via Luer-lock connection under laminar flow. The medium was transferred to the SmartBag bioprocess container, before it was reinstalled onto the SmartRocker platform. About 10 minutes prior the inoculation, single-use optical pH and DO Finesse TruFluor sensors were calibrated both by a one-point calibration.

The seeding procedure was realized in a laminar flow cabinet, where 1 L of the BY-2 inoculum suspension and 1 L fresh MSMOplus medium were transferred into the SmartBag bioprocess container by using an autoclaved funnel with male Luer-lock connector in order to achieve the desired cell concentration of 10 % pcv at the desired starting volume of 5 L. Subsequently, the SmartBag bioprocess container was reconnected to the control unit, the heating blanket for the exhaust air filter was installed, and all control loops for rocking motion, DO and temperature were started. The settings were defined as given below. About 60 min after the cell transfer, the first sample was taken using a specially designe sampling device (see Fig. 2).

Three individual runs with similar culture conditions (maximum rocking rates was lower in run#2 and run#3) were realized in order to evaluate the reproducibility of the results.

Culture conditions

The cultivation conditions were as follows:

Culture volume:	5 L
Rocking rate:	Run 1: 24 – 40 rpm Run 2&3: 24 – 36 rpm (increased step-wise daily after the sampling)
Rocking angle:	8 – 10 °
pH value:	non-regulated
Temperature:	26°C
Dissolved oxygen	50 %sat
Air flow rate:	0.25/0.1* slpm
Max. O ₂ flow rate:	0.25/0.1* slpm
Start cell concentration:	10 % pcv
Cultivation time:	10 days

* In run#1, flow rate was decreased to 0.1 slpm after the clogging of the outlet filter.

The bioreactor was agitated with increasing rocking rates and angles in order to ensure sufficient mixing of the culture broth, which became increasingly viscous (increase by a factor of 70 or greater [Eibl *et al.*, 2009])

because of biomass production and secretion of polysaccharides.

The dissolved oxygen concentration in the Finesse SmartRocker bioreactor was controlled solely by the addition of pure oxygen if the controller output of the DO regulation exceeded 70 %. No regulation of pH was realized. The temperature was controlled to 26°C via the integrated heater of the rocking platform.

Sampling and analysis

Samples were taken once a day with a specially designed sampling device mainly consisting of a Duran bottle (250 mL), which was sealed to the thermoplastic C-Flex tube (1/4 inch) of the SmartBag bioprocess container by using the Rewelder from ReedElectronics (see Fig. 2). Samples of about 40 - 100 mL were taken from the SmartBag bioprocess container via a vacuum pump. For in-process control (IPC), biomass growth was determined directly by analyzing fresh cell weight, dry cell weight and pcv. Indirect methods to determine the biomass growth included pH and conductivity measurements.




Fig. 2: Picture of the BY-2 tobacco cell culture in the SmartBag bioprocess container after seven days of culture. The specially designed sampling system is shown on the left side.

The pcv, which was defined as the volume ratio of packed biomass to the total suspension volume, was obtained by centrifugation of 2 x 10 mL cell suspension at 201 rcf for 15 min. For the determination of the fresh cell weight, 10 mL of the cell suspension was filtered for 3 min by using a filtration unit equipped with a round filter. The filtrated biomass was transferred to a dried (two days, 80°C) petri dish and the weight of the filtered biomass was calculated and converted to g·L⁻¹. The petri dish with the biomass was sealed by parafilm

and stored at -20°C. Subsequently, the petri dish containing the frozen biomass was placed in the drying oven for two days at 80°C. Conductivity and pH were determined in the supernatant generated during centrifugation for pcv determination using standard probes. For HPLC analysis of sucrose, glucose, fructose, ammonium and nitrate, 1.2 mL of cell free supernatant were frozen.

Results

During the experiments good agreement was obtained between the replicates (see Fig. 3). The highest deviation was found for the fresh dry weight, where all measurement points were within a range of ± 20 % around the mean value. However, very similar trends and excellent agreement of the final fresh weights (±4 %) were obtained. Therefore, only mean values of the offline data are given in the following.

The cells grew exponentially for three days with mean growth rates of 0.62±0.02 d⁻¹. A one-day lag-phase in run#1 was indicated by the fresh cell weight, but this was not confirmed by the packed cell volume (data not shown) and substrate consumption (see below). The determined growth rates are within typical ranges reported for wave-mixed and stirred single-use bioreactors, where growth rates between 0.54 d⁻¹ and 0.79 d⁻¹ have been reached [Raven *et al.*, 2010]. At day 4, the cells entered the stationary phase, where the biomass increased subsequently only slightly from 301±36 g·L⁻¹ to with 331±13 g·L⁻¹. Even though direct comparison is difficult because of different cultivation systems and operation conditions, the maximum biomass concentration (fresh weight 352±9 g·L⁻¹) is about 20 % lower than values reported by Kirchhoff (2012) for the Applikon® 3L stirred bioreactor (fresh weight 412 g·L⁻¹). Even higher fresh biomass concentrations of up to 470 g·L⁻¹ and 95 % pcv have been achieved in the Finesse SmartGlass™ 3L bioreactor (Finesse Solutions, Inc.) in our laboratory [Kaiser *et al.*, 2013]. Nevertheless, the cell viability remained high around 95 %, except of cultivation run#1, where the viability decreased to ≈ 90 % after 6 days of cultivation. This decrease may be explained by the fact that the exhaust filter was clogged. Due to the bag deformation, the cell suspension was less mixed. In subsequent cultivations, the exhaust air filter was changed at day 6, in order to prevent filter clogging.

As shown in Fig. 3b, sucrose was rapidly consumed in all cultivations, which led to an increase of glucose and fructose concentrations up to $11 \text{ g}\cdot\text{L}^{-1}$ (run#2, day 2) and $12 \text{ g}\cdot\text{L}^{-1}$ (run#3, day 1), respectively. Afterwards, the concentrations of both substrates decreased and were below the detection limit of $0.1 \text{ g}\cdot\text{L}^{-1}$ after four days of cultivation.

Typical profiles for the nitrate and ammonia ions were also found in all experiments. The nitrate and ammonium concentrations decreased during the exponential growth phase from the initial values of $2.0\pm 0.1 \text{ g}\cdot\text{L}^{-1}$ and $275\pm 40 \text{ mg}\cdot\text{L}^{-1}$ respectively. The minimum values of $17\pm 23 \text{ mg}\cdot\text{L}^{-1}$ (nitrate) and $3.67\pm 3.2 \text{ mg}\cdot\text{L}^{-1}$ (ammonium) were detected on day 5, before the concentrations of both ions increased, which may be attributed to cell lysis.

Due to the loss of ions in the culture broth, the conductivity decreased as well. The final value was $\approx 0.5\pm 0.1 \text{ mS}\cdot\text{cm}^{-1}$, whereas the initial conductivity was $5.1\pm 0.1 \text{ mS}\cdot\text{cm}^{-1}$. The offline determined pH values, increased over the entire cultivation from $\text{pH } 5.0\pm 0.4$ to $\text{pH } 7.5\pm 0.1$, since the pH value was non-regulated.

It should be emphasized that, based on our experiences using other single-use bioreactors with optical sensors, it is still challenging to accurately measure pH of plant cell cultures online. Often the pH value is below the measurement range, in particular at the beginning of the cultivation. Furthermore, some sensors are pre-calibrated by the manufacturer using a temperature of 37°C (instead of the typical 26°C for plant cell cultures). These factors may affect the reliability of the measurement. Indeed, no signal of the pH sensor was obtained during the first four days of the first experiment (run#1), which may be caused by the low pH of 5.4 (as determined offline) at the beginning of the cultivation, being at the lower limit of the sensor's measurement range [Finesse, 2014]. But, the sensors worked appropriately over the complete duration of the other two runs, even though recalibration was required almost daily (see Fig. 4a). Maximum deviation from offline to online pH value was 0.6 units. Recalibration was realized if the offline measured pH values of the samples differed more than 0.1 pH unit from the online data.

In contrast, the DO measurement worked well over almost the complete process. Some obstacles appeared for about 12 hours between the second and third day of cultivation, when the cells started to enter the

stationary phase. By the addition of pure oxygen in the gas supply and increase of rocking rate and angle, as shown in Fig. 4c, the DO remained among the setpoint and increased continuously after five days, so that an oxygen limitation can be excluded. Because of the surface aeration in the wave-mixed, no flotation was observed and no foam was produced (see Fig. 2), as often observed in sparged cultures [Raven *et al.*, 2010]. Hence, no addition of antifoam agent was required.

Nevertheless, rather high DO fluctuations between 10 %sat and 70 %sat were observed, which can be explained by the low oxygen mass transfer rates and probe response in the increasingly viscous culture broth. Perhaps, the DO controller settings should be adapted for further investigations.

After 5 days of cultivation, the dissolved oxygen started to increase because of the decreasing oxygen demands of the culture, which can be explained by the progressive cell death. At the end of the process, the culture broth was nearly oxygen-saturated.

In summary, it can be stated that *Nicotiana tabacum* cv. BY-2 suspension cells were successfully expanded using the novel wave-mixed Finesse SmartRocker bioreactor at 5 L scale. Experiments with other plant cell types in our laboratory are currently under preparation.

References:

J. Brändli, *et al.* (2012) Antikörperproduktion in Pflanzenzellen: Prozessentwicklung und -übertragung vom 50mL auf den 10L Masstab, *Biospektrum*, 2, 216-217.

R. Eibl, *et al.* (2009) Disposable bioreactors for plant liquid cultures at L-scale, *Engineering in Life Sciences*, 9, 156-164.

S.C. Kaiser, *et al.* (2013). Cultivation of *Nicotiana tabacum* cv. BY-2 suspension cells and IgG antibody production in the Finesse SmartGlass bioreactor. *Application Note*, ZHAW, unpublished.

J. Kirchoff (2012). Generation of highly productive polyclonal and monoclonal tobacco suspension lines from a heterogeneous transgenic BY-2 culture through flow cytometric sorting, *PhD-Thesis*. RWTH Aachen University.

N. Raven, *et al.* (2010). In: *Single-Use Technology in Biopharmaceutical Manufacture*, Eibl, R. & Eibl, D. (eds.), John Wiley & Sons, Inc., 251-261.

S. Werner, *et al.* (2010) Innovative, non-stirred bioreactors in scales from milliliters up to 1000 L for suspension cultures of cells using disposable bags and containers - a Swiss contribution, *Chimia*, 11, 819-823.

Finesse (2014). TruFluorPH specifications. Available online: www.finesse.com




Fig. 3: Profiles of offline data determined in a batch experiment with *Nicotiana tabacum* cv. BY-2 suspension cells. (a) Fresh cell weight and viability, (b) Concentrations of sugars sucrose, glucose and fructose, (c) Concentrations of nitrate and ammonium ions, (d) Conductivity and pH value.




Fig. 4: Profiles of online data determined in the batch experiment run#3 with *Nicotiana tabacum* cv. BY-2 suspension cells. (a) pH value measured online and offline – arrows indicated recalibration; (b) Dissolved oxygen and oxygen flow rate; (c) Rocking rate and angle. The arrows indicate recalibration of the pH measurement.