



Features

- Interfaces to Standard Ingold, B. Braun safety and sanitary ports
- FDA approved wetted materials
- Pg 13.5 threaded insertion for sensors
- Guarded and unguarded designs
- Material certs available

TruHousings™

TruHousings™ are designed to protect valuable electrodes and sensors when used in bioreactors or fermentors. These housings are manufactured from 316L S.S. and have been designed to interface with Standard Ingold, B. Braun safety and sanitary ports.

Specifications


Wetted Material	316L stainless steel, Viton o-rings
Sensor Connection	Pg 13.5
Port Connections	Standard, safety and 1½" sanitary
Pressure	6 bar maximum
Insertion	Guarded or unguarded

Housings (Side Mount, Guarded)

HSG-STD-025-G	fits 25 mm standard side port
HSG-SFT-025-G	fits 25 mm safety side port
HSG-SNY-075-G	fits 1½" sanitary side port

Housings (Side Mount, Unguarded)

HSG-STD-025-U	fits 25 mm standard side port
HSG-SFT-025-U	fits 25 mm safety side port



Guarded