

TruViu™ RDPD Part Numbers

I/O Tower

Family	Type	Tower ID	Gateway	Agitator	Sensors	Power Cord
TRV TruViu™	-XXXX	-XX	-XX	-YY	-P#D#X#X#	-AC
	-RDPD R&D, Process Dev	-IO I/O tower	-MB Modbus I/O gateway	-00 No Agitator	-P pH sensor, -D DO ₂ sensor	-US
			-FF Fieldbus I/O gateway	-AG Agitator Motor	A Auxilliary I/O	-UK
					-G Extra MFC, -L Extra pump -# 0 ... 1	-EU

Examples

TRV-RDPD-IO-MB-00-P1D1N0N0-US

TRV-RDPD-IO-FF-AG-P1D1G1A1-UK

Junction Box

Family	Type	Tower ID
TRV TruViu™	-XXXX	-XX
	-RDPD R&D, Process Dev	-JB Junction Box

2008-0317

Pump Tower



Family	Type	Tower ID	Pumps	Pumps
TRV TruViu™	-XXXX	-XX	-P#P#	-P#P#
	-RDPD R&D, Process Dev	-PM Pump Tower	-W3M# Watson Marlow 313 head, #	-W3M# Watson Marlow 313 head, #
		-SP Standalone Pump	# 0..4	# 0..4
			-WM323 Watson Marlow 323 pump	
			-WM520 Watson Marlow 520 pump	

Examples

- TRV-RDPD-PM-W3M3
- TRV-RDPD-SP-WM323