



Eat or be eaten.

Contaminants shouldn't have a chance.

Keep your media from being someone else's lunch. Our single-use sensors protect your bioprocess from predatory micro-organisms that can overrun cells and ruin your batch.

Finesse single-use sensors: SmartPuck pH+dO₂, TruFluor pH, TruFluor DO, and TruTorr block external contaminants. The sensors are inserted or welded into the bioreactor bag before sterilization. The readers never contact the media so the sterile barrier remains unbroken. TruFluor pH and TruFluor DO also have temperature-sensing capability, so no additional temperature sensor is needed.

All of the wetted materials are USP class VI and ACDF certified. Sensors are pre-calibrated, provide industry leading precision, and are plug-and-play in most any system.

Don't let other people's sensor claims leave a bad taste in your mouth.



Finesse Single-Use Sensors